



1st rapid exchange PTA balloon for intracranial stenosis



- Long soft distal tip with 0.016" thin profile
- Low balloon crossing profile: 0.023" - 0.026"
- Compatible with 0.014" guide wire
- "Zero profile" x-ray markers
- Total length: 160 cm
- Diameter sizes from 1.25 to 4.0 mm
Balloon lengths from 10 to 40 mm
- Nominal pressure at 6 ATM*
Rated burst pressure at 14 ATM*
- Non-/semi-compliant**


0.016" entry profile
embedded x-ray markers



*For pITA with diameter of 1.25 mm and 1.5 mm: Nominal pressure at 10 ATM and rated burst pressure at 18 ATM

**For pITA with diameter of 1.25 mm and 1.5 mm: Non-compliant; For all others: Semi-compliant

Rapid Exchange Neurovascular PTA Balloon Catheter



*The pITA with diameter of 1.25 mm and 1.5 mm have only one centrally located x-ray marker only.

REF	A: Balloon Diameter [mm]	B: Balloon Length [mm]	Nominal Pressure [atm]	Rated Burst Pressure [atm]
PITA-125-10	1.25	10	10	18
PITA-150-10	1.5	10	10	18
PITA-200-10	2.0	10	6	14
PITA-225-10	2.25	10	6	14
PITA-250-10	2.5	10	6	14
PITA-275-10	2.75	10	6	14
PITA-300-10	3.0	10	6	14
PITA-350-10	3.5	10	6	14
PITA-400-10	4.0	10	6	14
PITA-125-15	1.25	15	10	18
PITA-150-15	1.5	15	10	18
PITA-200-15	2.0	15	6	14
PITA-225-15	2.25	15	6	14
PITA-250-15	2.5	15	6	14
PITA-275-15	2.75	15	6	14
PITA-300-15	3.0	15	6	14
PITA-350-15	3.5	15	6	14
PITA-400-15	4.0	15	6	14
PITA-125-20	1.25	20	10	18
PITA-150-20	1.5	20	10	18
PITA-200-20	2.0	20	6	14
PITA-225-20	2.25	20	6	14
PITA-250-20	2.5	20	6	14
PITA-275-20	2.75	20	6	14
PITA-300-20	3.0	20	6	14
PITA-350-20	3.5	20	6	14
PITA-400-20	4.0	20	6	14

REF	A: Balloon Diameter [mm]	B: Balloon Length [mm]	Nominal Pressure [atm]	Rated Burst Pressure [atm]
PITA-200-25	2.0	25	6	14
PITA-225-25	2.25	25	6	14
PITA-250-25	2.5	25	6	14
PITA-275-25	2.75	25	6	14
PITA-300-25	3.0	25	6	14
PITA-350-25	3.5	25	6	14
PITA-400-25	4.0	25	6	14
PITA-200-30	2.0	30	6	14
PITA-225-30	2.25	30	6	14
PITA-250-30	2.5	30	6	14
PITA-275-30	2.75	30	6	14
PITA-300-30	3.0	30	6	14
PITA-350-30	3.5	30	6	14
PITA-400-30	4.0	30	6	14
PITA-200-35	2.0	35	6	14
PITA-225-35	2.25	35	6	14
PITA-250-35	2.5	35	6	14
PITA-275-35	2.75	35	6	14
PITA-300-35	3.0	35	6	14
PITA-350-35	3.5	35	6	14
PITA-400-35	4.0	35	6	14
PITA-200-40	2.0	40	6	14
PITA-225-40	2.25	40	6	14
PITA-250-40	2.5	40	6	14
PITA-275-40	2.75	40	6	14
PITA-300-40	3.0	40	6	14
PITA-350-40	3.5	40	6	14
PITA-400-40	4.0	40	6	14